

## Copyright

Shenzhen Comen Medical Instruments Co., Ltd.

Version: B00

No.: 046-00000248-00

Revision date: 11/2022

Product Name: Operating Table

Product Model: W5, W5A, W3 and W3A

Production Date: Refer to the label

Service Life: 10 years

## Statements

*All rights reserved by Shenzhen Comen Medical Instruments Co., Ltd.*

Shenzhen Comen Medical Instruments Co., Ltd. (hereinafter referred to as Comen or Comen Company) owns the copyright of this non-publicly published user manual and has the right to treat it as confidential information. This user manual contains proprietary materials protected by copyright law. No one shall disclose any content of the manual to others without the prior written consent of Comen. No part of the manual shall be photocopied, copied or translated to any other language.

This user manual does not imply the transfer of any proprietary right under patent law to any third party. Comen shall not be liable for any legal consequences arising out of breach of the patent law or infringement of any third-party rights.

This manual only provides reference information for operation, maintenance and repair of Comen products. Comen holds the final explanation right of this manual. Comen shall assume no liability for any errors contained herein or any incidental or indirect damages arising from provision, actual performance or use of this manual.

The contents contained in the manual may be changed without notice.

## Warranty

Comen will be responsible for the safety, reliability and performance of the product within the warranty period, if all of the following conditions are satisfied:

- The product is used in accordance with this manual.
- The product is installed, maintained or upgraded by Comen approved or authorized personnel.

- The storage and operating environments for the product shall comply with the recommended information and product specifications contained in this manual.
- The serial number label or manufacturing mark of the product is clearly legible.
- The damage is not caused by human factors.
- All replaceable components used for maintaining the device, accessories, and consumables are originally supplied by Comen or recognized by Comen.

Products conforming to Comen's warranty service policy can enjoy free services; for products not covered by the warranty scope, Comen will provide customer-paid services. If a product is transported to Comen for repair, the transportation expense (including customs fees) should be on the customer's account.

## **Return Procedure**

Please observe the following steps if it is necessary to return your product to Comen:

Obtain the return merchandise authorization (RMA): Contact the After-sales Service Department of Comen and inform Comen of the product serial number, which can be found on the product's nameplate. If the serial number is not clear and legible, the product return may not be accepted. Please indicate the product serial number and the date of manufacture and briefly describe the reason for return on the after-sales application form provided by the Department.

## **After-sales Service Department**

Name: After-sales Service Department of Shenzhen Comen Medical Instruments Co., Ltd.

Address: Floor 5 of Building 2, FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, Guangdong, 518106, P.R. China

Tel: +86-755-26431236 +86-755-86545386 +86-755-26074134

Fax: +86-755-26431232

Customer Service Tel: +86 4007009488

## **Preface**

This user manual provides details on the performance, operations and safety instructions relating to the W5, W5A, W3 and W3A operating tables. Please read it carefully and understand the content of this manual so as to ensure the safety of patients and operators.

This manual introduces the product according to the most complete configurations. The product you purchased may not have certain configurations or functions.

Please place the user manual near the operating table so that it can be easily and timely obtained when needed.

## **Intended Users**

This manual is suitable for professional clinical personnel or those who are expected to have knowledge and work experience in medical procedures, practices, and terminology necessary to monitor patients.

## **Illustrations**

All illustrations provided in this manual are for reference only. The menus, options, values and functions in the illustrations may not exactly match what you see in product.

BLANK

## Contents

<b>Chapter 1 Safety .....</b>	<b>1-1</b>
1.1 Safety Warning .....	1-1
1.2 Contraindications.....	1-6
1.3 Symbols .....	1-7
<b>Chapter 2 Product Overview .....</b>	<b>2-1</b>
2.1 Product Composition .....	2-1
2.2 Applicable To .....	2-1
2.3 Intended Use .....	2-1
2.4 Product Appearance .....	2-2
2.4.1 Front and Right Hand Views of W5 and W5A.....	2-2
2.4.2 Left Hand and Upward Views of W5 and W5A.....	2-3
2.4.3 Left Hand View of W3 and W3A.....	2-4
2.4.4 Hand Control .....	2-5
2.4.5 Side Control Panel .....	2-7
<b>Chapter 3 Installation and Disassembly .....</b>	<b>3-1</b>
3.1 Installation Requirements .....	3-1
3.1.1 Unpacking and Checking .....	3-1
3.1.2 Environmental Requirements.....	3-2
3.2 Installation of the Table .....	3-2
3.3 Installation and Disassembly of the Upper Back Plate .....	3-2
3.4 Installation and Disassembly of the Head Plate .....	3-3
3.5 Installation and Disassembly of the Leg Plate .....	3-4
3.6 Installation and Disassembly of the Mattress.....	3-5
3.7 Installation of Accessories .....	3-6
3.7.1 Clamp .....	3-6
3.7.1.1 Standard Circular Clamp .....	3-7
3.7.1.2 Standard Quadrate Clamp.....	3-7
3.7.2 Anesthesia Frame .....	3-8
3.7.3 Hand Support, Shoulder Support, Body Support, Foot plate, and Leg Support .....	3-9
3.7.4 Fixing Device for Tubes and Lines.....	3-10
3.8 Equipment Preparation .....	3-11
3.8.1 Connect the AC Power Supply .....	3-11
3.8.2 Protective Grounding .....	3-11
3.8.3 Condensation .....	3-12
3.9 Startup and Shutdown.....	3-12
3.9.1 Startup.....	3-12

3.9.2 Shutdown .....	3-13
<b>Chapter 4 Basic Operation .....</b>	<b>4-1</b>
4.1 Activate the Functions of the Operating Table .....	4-1
4.2 Longitudinal Adjustment Operation .....	4-1
4.3 Trendelenburg/ Reverse Trendelenburg Operation.....	4-1
4.4 Lateral Lift Operation.....	4-2
4.5 Operation of the Head Plate.....	4-2
4.6 Operation of the Upper Back Plate.....	4-3
4.7 Operation of the Back Plate.....	4-4
4.8 Operation of the Leg Plate.....	4-4
4.9 Lifting Operation.....	4-5
4.10 Forward/Backward Buckling Operation with One Key .....	4-5
4.11 Operation of the Waist Elevator .....	4-6
<b>Chapter 5 Fault Handling .....</b>	<b>5-1</b>
<b>Chapter 6 Battery .....</b>	<b>6-1</b>
6.1 Overview .....	6-1
6.2 Battery Installation .....	6-2
6.3 Optimization and Inspection of the Battery Performance.....	6-2
6.4 Battery Recycling .....	6-3
<b>Chapter 7 Cleaning and Disinfection .....</b>	<b>7-1</b>
7.1 Overview .....	7-1
7.2 Cleaning and Disinfection of the Operating Table .....	7-2
7.3 Cleaning and Disinfection of Accessories .....	7-3
<b>Chapter 8 Maintenance .....</b>	<b>8-1</b>
8.1 Maintenance inspection.....	8-1
8.2 Maintenance Plan.....	8-2
8.2.1 Maintenance .....	8-2
8.2.2 Service Life of Accessories.....	8-2
<b>Appendix I Accessories.....</b>	<b>I-1</b>
<b>Appendix II Product Specification .....</b>	<b>II-1</b>
<b>Appendix III EMC .....</b>	<b>III-1</b>
<b>Appendix IV Toxic and Hazardous Substances or Elements in the Product .....</b>	<b>IV-1</b>
<b>Appendix V Consideration for Environmentally Conscious Design.....</b>	<b>IV-1</b>

### 1.1 Safety Warning

---



#### Warning

- It indicates serious consequences, adverse events or situations that endanger safety. Failure to follow the warning will result in serious personal injury or death to the user or patient.



#### Caution

- It indicates a potentially dangerous or unsafe operation that, if not avoided, may cause minor personal injury, product failure or damage or property, or cause more serious injury in the future.



#### Attention

- It emphasizes important precautions and provides instructions or explanations for better use of the product.



#### Warning

- Be sure to read and understand the entire manual before using the operating table. As with all medical devices, attempts to use this equipment without a thorough understanding of its operation can result in injury to the patient or user.
- The operating table can only be operated by trained personnel. Their training records must be archived. Any unauthorized or untrained personnel are not allowed to perform any operations.
- Be sure to keep this manual near the operating table for convenience.
- Pay attention to various warnings or prompts on the surface of the equipment to ensure the safety of operators and the normal operation of the equipment.
- This product was delivered without disinfection. It must be cleaned and disinfected before the first use.

- This equipment can only be used for the specified intended use and only allowed to be used by trained and qualified personnel.
- Before use, the user must check whether the equipment and its accessories can work normally and safely.
- This equipment can only move on flat ground and cannot be pushed to pass over steps or thresholds to avoid the risk of tipping over.
- This equipment cannot be used in the presence of flammable anesthetic gas mixed with air, oxygen or nitrous oxide.
- Please carefully plug in and route the power cord and various accessory cables to avoid the patient being entangled or suffocated, cable entanglement, or electrical interference.
- Do not place the power plug that disconnects the equipment from the AC power supply in a location that is difficult for the operator to access.
- During use the operator shall always observe whether the equipment is normally connected to the power supply so as to avoid adverse events such as shutdown caused by loose power supply connection.
- Be sure to connect the protective ground terminal.
- This equipment can only be connected to a power socket with protective grounding. If the power socket is not connected to a ground conductor do not use the socket and use a rechargeable battery for power supply instead.
- If there is any doubt about the installation of the external protective ground terminal or the integrity of the wiring, the internal battery must be used to run the operating table.
- If there is an abnormality in the power supply of the equipment or during the use of the equipment, please use the emergency stop switch to terminate the current abnormal state.
- The equipment interconnected with this device should form an equipotential body (effectively connected to the protective ground).
- When using the operating table, be sure to cover the mattress with a sterile sheet to avoid infection. Be sure to avoid using wet sheets.
- Make sure to use this equipment on an anti-static floor.
- Please use the anti-static mattress provided by Comen. If other mattresses are used, it may cause injury to the patient or damage to the equipment.
- Do not use alcohol-containing detergents/disinfectants where high-frequency surgical equipment may be used to avoid fire.



- When using equipment such as a high-frequency electrotome, defibrillator, etc., if the patient touches the metal parts of the operating table (including accessories), or lies on wet sheets or a conductive mattress, there may be a risk of burns. During the operation, prevent the patient from touching the metal parts of the operating table (including accessories), and ensure that dry sheets and non-conductive mattresses are used.
- Electromagnetic fields can affect the performance of this equipment. Therefore, other equipment used in the vicinity of this equipment must meet the corresponding EMC requirements. For example, mobile phones, walkie-talkies, MRI and X-rays may be sources of interference as they all emit high-intensity electromagnetic radiation.
- The mass equivalent filtration of the operating table shall not exceed 1mmAL.
- When moving the operating table please make sure that no objects are placed on the table surface and ensure that the passage is free of obstacles to avoid damage or collision.
- When moving the operating table please use both hands to ensure the direction of movement to avoid collisions.
- When adjusting the table surface to fall down or tilt or folding down the head/foot plate, please be careful that the table surface or head/foot plate does not collide with the base or the ground.
- When adjusting the operating table please be mindful of hanging parts, sheets, etc. colliding or getting stuck.
- Before unlocking the operating table be sure to grasp the operating table firmly to prevent unexpected movement of the operating table after unlocking.
- Please pay attention to the location of pinch points to avoid injury to the patient or operator.
- Please place or transfer patients under the guidance of medical personnel.
- Before placing the patient make sure that the width of the operating table is appropriate, otherwise it may cause personal injury.
- Before placing the patient make sure that the casters of the operating table are locked to prevent unexpected movement.
- Do not move the operating table when the patient is on the table.
- When using the operating table be sure to confirm that it is within the maximum load-bearing range (normal position or interchange position of the head and leg plates).
- If accessories are to be used during the operation, make sure that the sum of the weight of the accessories and the patient does not exceed the maximum load-bearing range (normal position or interchange position of the head and leg plates).
- When the operating table is in some special positions there are certain restrictions on the load

bearing. Please be sure to consider the maximum load bearing and related restriction conditions when using the table. If in doubt, please consult the manufacturer or local agent before use.

- Before transferring the patient to the operating table make sure that the operating table is locked and the entire table surface is in the best balanced position: that is, all the table surfaces are level, and the entire table surface is at the extreme position of longitudinal slide (foot). Do not transfer the patient in other positions, otherwise the table body may tip over.
- When transferring the patient please transfer from the side. Improper operation may cause the operating table to tip over.
- Improper placement of the patient will cause the operating table to tip over. When placing the patient make sure that the direction is correct. Do not place the patient's body on the leg plate.
- When placing the patient make sure that the patient's center of gravity is as close as possible to the center of the column, otherwise there is a danger of tipping.
- Incorrect posture may cause damage to the patient's physical functions, especially when the patient is in a lateral tilt or trendelenburg / reverse trendelenburg position at a large angle. When using the operating table make sure that the patient's position is correct, and pay attention to check whether the position is offset at any time.
- While ensuring the correct position attention should be paid to the patient's condition to avoid damage to the patient's respiratory system, nervous system and circulatory system.
- When using the operating table, be sure to select the correct operation according to the patient's position, otherwise it will cause personal injury.
- The weight of different patients will affect the movement speed of the table surface, and thus please be careful during the adjustment process.
- When the patient weighs above 150kg, the trendelenburg / reverse trendelenburg angle shall not exceed 10°, and the lateral tilt angle shall not exceed 5°, so as to prevent the patient from slipping and causing personal injury.
- When the table surface is tilted toward the foot please be careful of the collision between the head plate and the patient's head when folding the back plate.
- Before using the longitudinal slide function make sure that the surface of the operating table remains level.
- After using the longitudinal slide function you can perform lateral tilt and trendelenburg / reverse trendelenburg, but not the fold up/ fold down and other movements of the back plate. Please always pay attention to the patient's position during operation to avoid collision.
- When adjusting the position of the operating table or patient, especially when the patient is in a

lateral tilt or trendelenburg / reverse trendelenburg position at a large angle, make sure that the patient has been properly fixed, otherwise it may cause personal injury or equipment damage.

- If you use the longitudinal slide function to perform X-ray radiography on the patient before a surgical operation, after radiography, be sure to restore the table to a balanced position (translate the foot to the extreme position) before performing the operation.
- When adjusting the head plate, back plate, leg plate and other parts of the operating table, be sure to be careful of the joints to avoid pinching the patient or operators.
- When adjusting or moving the operating table be sure to pay attention to the position of the operating table, accessories, patients and other items to avoid personal injury or equipment damage due to collisions.
- When the table surface is being adjusted to a lower position, while operating the table surface, please be careful of the collision or pinching of the table surface with the base or the ground.
- When the head/leg plate is folded down, please be careful of the collision or pinching of the head/leg plate with the base or the ground.
- When tilting the table surface or adjusting the leg plate make sure that the table surface does not collide with the column or base.
- When adjusting the surface of the operating table make sure that the wire controller cable is not trapped by the table surface joint.
- When adjusting the operating table please be careful that the hanging parts (such as the leg plate, hand support, etc.) do not collide with objects in the operating room.
- When adjusting the operating table please be careful to prevent accessories, sheets, etc. from getting stuck in the moving parts of the operating table to avoid equipment damage.
- When the operating table is not in use the user must ensure that the equipment is locked and the entire table surface is in the optimal position: that is, all the table surfaces are level, the entire table surface is at the extreme position of longitudinal slide (foot), and the table surface is lowered to the extreme position.
- The safety and effectiveness of the operating table can only be guaranteed if it is properly maintained. Be sure to maintain the operating table according to the manufacturer's instructions.
- Do not disassemble the housing of the equipment to avoid possible electric shock. Any maintenance and upgrade of the equipment must be carried out by service personnel trained and authorized by the Comen company.
- Packaging materials must be disposed of in accordance with local laws and regulations or hospital waste disposal rules and regulations. Packaging materials must be placed out of the reach of children.

- The operating table can only be installed by authorized personnel of the Comen company.
- The operating table must be thoroughly inspected for safety by professionals every year.
- In order to avoid environmental pollution please dispose of discarded equipment, old batteries or accessories in accordance with relevant local regulations or hospital requirements.



#### **Caution**

- In order to avoid equipment damage and ensure patient safety please use the accessories specified in this manual.
- Please install or move the equipment properly to prevent the equipment from falling, colliding, being damaged by strong vibration or other external mechanical forces.
- Before the equipment is powered on please confirm whether the mains power supply meets the requirements of the equipment nameplate label or the supply voltage and frequency specified in the User manual.
- In order to protect the environment disposable accessories must be recycled or properly disposed of.



#### **Attention**
























- Please install the equipment in a place convenient for observation, operation and maintenance.
- This manual introduces the product according to the most complete configurations. The product you purchased may not have certain configurations or functions.
- Please place the user manual near the equipment so that it can be quickly and easily obtained when needed.
- This equipment cannot be used at home.
- The equipment can only be used by one patient at a time.
- In normal use, the operator's position shall be within one meter of the equipment.

## **1.2 Contraindications**

None.





## 1.3 Symbols

### ● Instrument symbols

Symbol	Description	Symbol	Description
	Caution!		Refer to Instructions for Use
	Complies with Medical Device Regulation (EU) 2017/745		Authorized representative in the European Community
	Type B applied part		Protection against harmful effects of splashing water (overall unit)
	Protective grounding		Non-ionizing radiation
	Anti-pinch		Warning!
	Alternating current		Mains power input (national standard)
	Battery symbol		RS232 interface
	Serial number		Date of Manufacture
	Environmental protection		Classified collection of electrical and electronic equipment
	The environmental protection service life of electronic products is 20 years		Manufacturer
	Temperature limits		Humidity limitation
	Atmospheric pressure limitation		

Description: Please refer to Chapter 2 for the button symbols and functions of the operating table's controller panel.

### ● Packaging symbols

Symbol	Description	Symbol	Description
	Do not stack		Do not roll
	This way up		Fragile, handle with care

	Keep away from rain		Center of gravity
---	---------------------	---	-------------------



### Attention

- Please refer to the actual equipment and the label on the packaging for the equipment symbol.

---

## Chapter 2 Product Overview

---

The design of the equipment conforms to the relevant Chinese and international safety standards of medical electrical equipment.

### 2.1 Product Composition

---

The operating table is mainly composed of the table (including the support component, transmission component and electric control component) and accessories (optional). Among which, W5 and W5A table surfaces can be divided into five segments: head plate, upper back plate, back plate, hip plate and leg plate; W3 and W3A table surfaces can be divided into four segments: head plate, back plate, hip plate and leg plate.

The operating table can be equipped with general operating table accessories such as an anesthesia frame, hand support, leg support, foot plate, body support, shoulder support, clamp, prone position pad, body strap, and wire manager.

### 2.2 Applicable To

---

Patients in hospitals.


### 2.3 Intended Use

---

For clinical operations.

## 2.4 Product Appearance


### 2.4.1 Front and Right Hand Views of W5 and W5A






No.	Name	No.	Name
(1)	Head plate	(2)	Upper back plate
(3)	Back plate	(4)	Hip plate
(5)	Leg plate	(6)	Side rail
(7)	Fixed handle for leg plate	(8)	Horizontal locking wrench for leg plate
(9)	Column	(10)	Caster
(11)	Base	(12)	Fuse
(13)	AC power supply interface	(14)	Protective ground terminal
(15)	RS232 interface for wall screen control	(16)	Emergency stop switch

### 2.4.2 Left Hand and Upward Views of W5 and W5A



(1)	Hand Control	(4)	Adjustment lever for leg plate
(2)	Adjustment lever for head plate	(5)	Side control panel
(3)	Hand Control interface	(6)	Foot control interface

### 2.4.3 Left Hand View of W3 and W3A




(1)	Head plate	(2)	Back plate	(3)	Waist Elevator
(4)	Hip plate	(5)	Leg plate	(6)	Waist elevator rocker



#### Warning

- All analog and digital equipment connected with this device must be products certified by the specified IEC standards (such as IEC 60950 IT equipment standard and IEC 60601-1 Medical device standard). Moreover, all configurations shall follow the valid version of the IEC 60601-1 system standard. The personnel connecting additional equipment to the input/output signal ports are responsible to ensure that they configure the entire system appropriately so that it meets the requirements of the IEC 60601-1 standard. If there are any doubts, please contact our company.
- When signal interfaces such as the serial port are connected to multiple pieces of equipment at the same time, the total leakage current caused shall comply with IEC 60601-1.

## 2.4.4 Hand Control



No.	Name	Description
(1)	AC power indicator	<ul style="list-style-type: none"> <li>● Steady white light: The operating table is connected to an AC power supply.</li> <li>● Light off: The operating table is not connected to an AC power supply.</li> </ul>
(2)	Battery indicator	<ul style="list-style-type: none"> <li>● White light flashes: The operating table is powered by battery.</li> <li>● Yellow light flashes: The battery is low.</li> <li>● Steady white light: The battery is being charged.</li> <li>● Light off: The battery is fully charged.</li> </ul>
(3)	Locked indicator	<ul style="list-style-type: none"> <li>● Steady white light: The operating table is locked and cannot be moved.</li> <li>● Light off: The operating table is unlocked and can be moved.</li> </ul>

(4)	Unlocking indicator	<ul style="list-style-type: none"> <li>● Steady yellow light: The operating table is unlocked and can be moved.</li> <li>● Light off: The operating table is locked and cannot be moved.</li> </ul>
(5)	Shutdown	
(6)	Startup	The white indicator light is on when the machine is turned on, and the indicator light is off when it is turned off (the power will be turned off after 60s~70s without operation)
(7)	Lock system	Press the button for 3 seconds to lock the equipment.
(8)	Unlock system	Press the button for 7 seconds to unlock the equipment.
(9)	Longitudinal Slide (foot)	
(10)	Longitudinal Slide (head)	
(11)	Rise	
(12)	Fall	
(13)	Reverse Trendelenburg	
(14)	Trendelenburg	
(15)	Lateral tilt (right)	
(16)	Lateral tilt (left)	
(17)	Back plate fold up	
(18)	Back plate fold down	
(19)	Forward buckling	
(20)	Backward buckling	
(21)	"0" button	Press it to return to the original position with one push of the "0" button.

**Warning**

- When the patient's position adjustment is complete, be sure to turn off the hand control to avoid mis-operation.


**Attention**

- When the hand control is idle for 60 ~ 70s, it will turn off automatically.
- Before each use of the operating table check and confirm that all functions of the operating table can

operate normally. Operate all the function buttons one by one, each operation lasting more than 3 seconds, and finally press the “0” button to adjust the table surface to the original position.

- For the hand control of the W3 and W3A operating tables, longitudinal slide (head), longitudinal slide (foot), forward/backward bucking and “0” button are optional functions.

#### 2.4.5 Side Control Panel



No.	Name	Description
(1)	Locking indicator	<ul style="list-style-type: none"> <li>● Steady white light: The operating table is locked and cannot be moved.</li> <li>● Light off: The operating table is unlocked and can be moved.</li> </ul>
(2)	Unlocking indicator	<ul style="list-style-type: none"> <li>● Steady yellow light: The operating table is unlocked and can be moved.</li> <li>● Light off: The operating table is locked and cannot be moved.</li> </ul>
(3)	Lock system	Press the button for 3 seconds to lock the equipment.
(4)	Unlock system	Press the button for 7 seconds to unlock the equipment.
(5)	Longitudinal slide (head)	
(6)	Longitudinal slide (foot)	
(7)	Fall	
(8)	Rise	
(9)	Trendelenburg	

(10)	Reverse Trendelenburg	
(11)	Lateral tilt (left)	
(12)	Lateral tilt (right)	
(13)	Back plate fold down	
(14)	Back plate fold up	
(15)	Backward buckling	
(16)	Forward buckling	
(17)	"0" button	Press to return to the original position with one push of the "0" button.
(18)	Startup	The white indicator light is on when the machine is turned on, and the indicator light is off when it is turned off (the power will be turned off after 60s~70s without operation)
(19)	Shutdown	
(20)	Battery indicator	<ul style="list-style-type: none"> <li>● White light flashes: The operating table is powered by battery.</li> <li>● Yellow light flashes: The battery is low.</li> <li>● Steady white light: The battery is being charged.</li> <li>● Light off: The battery is fully charged.</li> </ul>
(21)	AC power indicator	<ul style="list-style-type: none"> <li>● Steady white light: The operating table is connected to an AC power supply.</li> <li>● Light off: The operating table is not connected to an AC power supply.</li> </ul>

---

## Chapter 3 Installation and Disassembly

---



### Attention

- In order to ensure the normal operation of the operating table please read this Chapter and the safety chapter before use, and install it as required.
- This product was delivered without disinfection. It must be cleaned and disinfected before the first use.



### Warning

- Be sure to check whether each part is in good condition, and replace it if it is damaged.
- When disassembling/installing parts of the operating table, be sure to hold the parts using both hands firmly to avoid equipment damage or personal injury caused by slipping.
- Be careful of pinching your hand when disassembling/installing each part.
- Do not install the hand support, anesthesia frame, and infusion frame on the side rail of head plate, so as to avoid unexpected movement of the head plate and thus cause equipment damage or personal injury.
- If you need to use accessories during the operation you must ensure that the sum of the weight of the accessories installed on each part and the weight of the patient's limbs does not exceed the maximum load bearing of the part.
- Be sure to use the specified mattress to avoid injury to the patient during the movement of the table surface.

## 3.1 Installation Requirements

---

### 3.1.1 Unpacking and Checking

---

Carefully take the device and its accessories out of the packing box, and save the packaging materials for future transportation or storage. Please count the accessories according to the packing list. Check for any mechanical damage. Check all external wiring and cables, and some inserted accessories. In case of any

problems, please contact our After-sales Department or agent immediately.



### **Warning**

- **If any damage is found, please contact the relevant personnel of the hospital or the after-sales Service Department of Comen.**

### **3.1.2 Environmental Requirements**

The working environment of the operating table shall comply with the requirements specified in this manual.

The ambient temperature beyond the specified ranges may affect the accuracy of the equipment and cause damage to the components and wiring.

The operating table shall be used in an environment reasonably free from vibration, dust, corrosive or explosive gases, extreme temperatures, humidity, etc.

## **3.2 Installation of the Table**

- a) Place the table in a suitable position;
- b) Press and hold the **[System Lock]** button on the side control panel to lock the table on the ground.
- c) If the table is unstable, please check whether the ground is flat and adjust the position of the table to ensure that the table is firmly locked to the ground.
- d) Insert the connecting plug of the hand control into the socket located on the column, and then hang the hand control on the side rail.


## **3.3 Installation and Disassembly of the Upper Back Plate**

The back plate is inserted into the table, and no installation is required. Just install the upper back plate.

- 1) The installation steps of the upper back plate are as follows:
  - a) As shown in the figure below, hold the upper back plate firmly with both hands, and align the quick release pins (3) with the quick release slots (2) on the left and right sides of the back plate;
  - b) Push the upper back plate into the slot firmly, the quick release button (1) of the back plate pops up, and you hear a “click” sound;
  - c) Tighten the knob (4);



- d) Check whether the installation is safe and reliable.




## 2) Disassembling the upper back plate:

- a) Before disassembling the back plate, the head plate must be disassembled first;
- b) Hold the upper back plate firmly with both hands, and press the quick release button (1) of the upper back board on the back plate; push the upper back plate in the direction of the back plate firmly with both hands, and the quick release button of the head plate on the back plate pops up;
- c) Pull out the upper back plate forcefully.

## 3.4 Installation and Disassembly of the Head Plate

### 1) The installation steps of the head plate are as follows:

- a) As shown in the figure below, hold the head plate firmly with both hands, and align the quick release pins (3) with the quick release slots (1) on the left and right sides of the upper back plate;
- b) Push the head plate into the slot firmly, the quick release button (2) of the head plate on the upper back plate pops up, and you hear a “click” sound;
- c) Check whether the installation is safe and reliable.




2) Disassembling the head plate:

- a) Hold the head plate firmly with both hands, and press the quick release buttons (2) of the head plate on the left and right sides of the upper back plate with your fingers;
- b) Pull out the head plate forcefully.

### 3.5 Installation and Disassembly of the Leg Plate

1) The installation steps of the leg plate are as follows:

- a) As shown in the figure, hold the joints of the leg plate with both hands, and align the quick release pin (3) with the quick release slot (2) of the hip plate;
- b) Push the plate forcefully into the slot, lift the fixed handle (1) for the leg plate and turn it tightly;
- c) Check whether the installation is safe and reliable.



2) Disassembling the leg plate:

- a) Hold the joints of the leg plate and loosen the fixed handle (1) for the leg plate with the other hand;
- b) Pull out the leg plate.

### 3.6 Installation and Disassembly of the Mattress



#### Warning

- Be sure to check whether the mattress is aging or damaged before use. If there is any damage, replace it immediately.
- Do not use a mattress with a damaged surface, otherwise liquid penetration will occur, which does not meet the hygiene requirements.
- When the patient is lying on the mattress for disinfection do not allow the disinfectant to accumulate under the patient's body.
- When installing a memory sponge mattress make sure that the Velcro on the mattress and the one on the table surface are in complete alignment, otherwise the mattress may slip off and cause personal injury.
- When you install the mattress, damaged, loose or wet Velcro will not guarantee the safety and reliability of the mattress installation, therefore please replace it immediately.
- When using the mattress be sure to cover the mattress with a sterile sheet to avoid infection.




#### Attention

- The mattress shall be placed horizontally.
- Before installing the mattress make sure that the table surface is dry.

This operating table is equipped with a memory sponge mattress, which is fixed on the surface of the operating table by Velcro. Each mattress can be disassembled.

- 1) The installation steps of the mattress are as follows:
  - a) Place the mattress on the table surface and make sure that the Velcro strips under the mattress are completely attached to those on the table surface, and press them firmly together.
  - b) Check whether the mattress installation is safe and reliable.



2) Disassembling the mattress:

Grasp one end of the mattress with both hands and pull up the mattress firmly.

### 3.7 Installation of Accessories

---

#### **Warning**

- The accessories are part of the operating table system which can only be used in conjunction with the operating table and accessories of the Comen company. Do not apply the accessories to other manufacturers' products, otherwise this may cause personal injury or equipment damage.
- Please check the condition of the accessories carefully before each use. If you find signs of aging or damage be sure to stop using the accessories and notify the relevant personnel to check and repair the accessories, and replace them if in bad condition.

#### 3.7.1 Clamp


---

The clamp is used to fix the operating table accessories on the side rail of the operating table. The operating table can be equipped with two types of clamps: standard circular clamp and light quadrangle clamp, which are used to install accessories equipped with round and square support rods. When using a clamp, be sure to use the suitable one.

#### **Warning**

- Only the recommended clamp can be applied to the side rail of the table.
- When installing the clamp, make sure that the clamp body is completely fixed on the side rail of the operating table without any tilt, otherwise the accessories may slip or fall during use, which may cause personal injury or equipment damage.
- After installation or adjustment be sure to firmly tighten the puller handle and clamp handle, otherwise the accessories may slip or fall during use, which may cause personal injury or equipment damage.

#### 3.7.1.1 Standard Circular Clamp




- (1) Clamp body
- (2) Clamp hole
- (3) Puller handle

Installation and usage methods:

- a) First unscrew the puller handle (3), and clamp the clamp body (1) on the side rail of the operating table;
- b) Insert the accessory to be fixed into the clamp hole (2);
- c) Tighten the puller handle to fix the accessory to the clamp hole and to lock the clamp body on the side rail;
- d) Check whether the installation is safe and reliable.

#### 3.7.1.2 Standard Quadrate Clamp




- (1) Clamp hole
- (2) Clamp body
- (3) Clamp handle

Installation and usage methods:

- a) Clamp the clamp body (2) on the side rail of the operating table;
- b) Insert the accessory to be fixed into the clamp hole (1), ensuring that the accessories are inside the side rail;
- c) Tighten the clamp handle (3), and lock the accessory on the clamp.
- d) Check whether the installation is safe and reliable.

### 3.7.2 Anesthesia Frame

The operating table can be equipped with a standard anesthesia frame for anesthesia curtains.



Installation and usage methods:

- a) Choose a suitable position and clamp the body of the circular clamp on the side rail of the operating table;
- b) Turn the tightening handle of the clamp to lock the clamp body on the side rail;
- c) Insert the anesthesia frame rod (1) into the clamp hole, adjust the rod up and down to adjust the anesthesia frame to a suitable height;
- d) Tighten the clamp handle to lock the anesthesia frame on the clamp;
- e) Hang the anesthesia curtain on the anesthesia frame;
- f) Check whether the installation is safe and reliable.




**Warning**

- The rod of this anesthesia frame is round and can only be used in conjunction with a circular clamp. Do not install this anesthesia frame with a quadrate clamp, otherwise this may cause personal injury or equipment damage.
- When installing the anesthesia frame, make sure that the anesthesia frame rod is vertically inserted into the clamp without any tilt, otherwise the anesthesia frame may slip or fall during use, which may cause personal injury or equipment damage.

### 3.7.3 Hand Support, Shoulder Support, Body Support, Foot plate, and Leg Support

The operating table can be equipped with rotatable hand support, shoulder support, body support, foot plate, and leg support, for placing and fixing the patient's body parts, one on each side.

The composition structure of the rotatable hand support is shown in the figure below:



Installation and usage methods:

- Choose a suitable position and clamp the clamp on the side rail of the operating table;
- Pull down the wrench (2) to lock the clamp on the side rail;
- Pull up the handle (5) to adjust the hand support to a suitable height;
- Fix the patient's arm on the hand support with a body strap (4);
- Check whether the installation is safe and reliable.



#### Warning

- The standard leg support can only be used in conjunction with the circular clamp. The standard shoulder support, standard body support, and standard foot plate can only be used in conjunction with the quadrate clamp. Do not fit any support into an incorrect clamp, otherwise it may cause

personal injury or equipment damage.


- When installing make sure that the part is inserted vertically in the clamp without any tilt, otherwise it may slip or fall during use, which may cause personal injury or equipment damage.
- Do not rest on the parts during the operation, otherwise it will cause permanent injury to the patient's arm. During the operation, be sure to check whether the patient's body is in the correct position.
- Be sure to adjust the hand support to be parallel to the table surface to prevent the patient from being injured when moving the patient.
- When using the hand support, shoulder support, body support, leg support and foot support, be sure to cover it with sterile sheets to prevent infection.



#### Attention

- The installation method of each accessory is similar, so you can refer to the hand support installation method when installing other accessories.

### 3.7.4 Fixing Device for Tubes and Lines



(1)	Clamp	(2)	Angle Adjustment Knob 1	(3)	Angle Adjustment Knob 2
(4)	Tube Clamp 1	(5)	Tube Clamp 2	(6)	Line Clamp

Installation and usage methods:



- a) Choose a suitable position and clamp the circular clamp on the side rail of the operating table;
- b) Turn the puller handle of the clamp to lock the clamp on the side rail;
- c) Insert the support rod of the fixing device into the clamp hole, and adjust the rod up and down to find a suitable height;
- d) Tighten the clamp handle, and lock the hand support on the clamp;
- e) Insert the cables into the line clamp;
- f) Check whether the installation is safe and reliable.

## 3.8 Equipment Preparation

### 3.8.1 Connect the AC Power Supply

Steps to connect the AC power supply:

- 1) Make sure that the AC power supply meets the specifications specified in this manual;
- 2) Use the power cord that was supplied with the operating table to connect one end of the power cord to the power connector on the equipment and the other end to a grounded electrical outlet;



#### Attention

- **Connect the power cord to the grounded hospital outlet of the hospital.**
- **When there is a battery configuration the battery must be charged after the equipment is transported or stored. If the equipment is turned on without connecting to AC power supply the equipment may not work normally due to insufficient battery power. After the AC power supply is connected the battery can be charged regardless of whether the equipment is turned on or not.**

### 3.8.2 Protective Grounding

In order to protect patients and operators the metal housing of the operating table must be grounded. Therefore, the operating table is equipped with a three-wire power cord. After plugging it into a matching three-wire socket, ground the equipment through the ground wire (protective grounding) in the power cord. If there is no three-wire socket, consult the electrical management personnel of the hospital.



#### Warning

- **Do not connect the three-wire power plug of this equipment to a two-wire socket.**
- **Be sure to connect the protective ground terminal.**
- **Make sure that the housing of the operating table is safely and reliably grounded.**

The equipment provides a protective ground terminal. Before using the equipment, be sure to connect the protective ground terminal. If it is unclear whether a particular combination of devices is dangerous (for example, danger from the accumulation of leakage current), the user shall consult the relevant manufacturer or other experts in this field to ensure that the device is used safely and will not be damaged by the recommended combination.

### 3.8.3 Condensation

---

When using this device please ensure that the equipment is free of condensation. When equipment is transferred from one room to another, condensation may form. This is because the equipment is exposed to humid air and different temperatures. In order to avoid potential risks (for example, condensed liquid can cause a short circuit on the circuit board), leave the instrument to dry before use.

Note: Condensation is the condensation of a gas or liquid when it cools. For example, steam turns to water when it cools, and water turns to ice when it cools. The lower the temperature, the faster the condensation forms.

## 3.9 Startup and Shutdown

---

### 3.9.1 Startup

---

This equipment can be powered by the built-in battery or an AC power supply.

- 1) Before startup check whether all parts of the operating table are connected correctly and whether there is any mechanical damage;
- 2) Check whether the operating table can be turned on normally:
  - Connect the power cord to the AC power supply, turn on the emergency stop switch (it automatically pops up when rotated), and after the power button is pressed the white AC power indicator will light up;
  - Press the ON button on the controller panel to start the equipment, and the indicators of the controller will illuminate.

**Attention**

- Check all available functions to make sure that the operating table is functioning properly.
- The battery of this equipment is optional. If a battery is fitted to this equipment, please charge the battery after each use to ensure that there is sufficient power reserve.
- It is recommended that the time interval of startup and shutdown is greater than 1 minute to prevent shortening the service life of the equipment.

### 3.9.2 Shutdown

---

If you need to turn off the operating table, please follow the steps below:

- 1) Confirm that the adjustment of the operating table is completed, and press the OFF button on the controller panel;
- 2) The hand control or side control panel will automatically shut down after 60~70s without operation;
- 3) After confirming that you can stop using the operating table, turn off the power switch to turn off the AC power indicator on the controller, and disconnect the AC power cord.
- 4) When the operating table will not be used for long periods, please press the emergency stop switch.

### 4.1 Activate the Functions of the Operating Table

---

All functions of the operating table can be activated through the hand control or the side control panel, and some parts can also be adjusted manually.

Operating steps of the hand control or the side control panel:


- 1) Turn on the emergency stop switch (automatically pops up when rotated);
- 2) Press the ON button to start the operating table;
- 3) Press the Startup button and the corresponding function button at the same time and release it when you reach the intended position, and the table will automatically stop.

### 4.2 Longitudinal Adjustment Operation

---

Press the Startup button and the Longitudinal Adjustment function button at the same time to adjust the length of the operating table.

The adjustment range is shown in the following figure.




### 4.3 Trendelenburg/ Reverse Trendelenburg Operation

---

Press the Startup button and the Trendelenburg/ Reverse Trendelenburg function button at the same time to adjust the operating table to stay in either position.


The adjustment angle range is shown in the following figure.



## 4.4 Lateral Lift Operation


Press the Startup button and the Lateral Lift function button at the same time to make the operating table lean to the left or right.

The adjustment angle range is shown in the following figure.




## 4.5 Operation of the Head Plate

- Firmly hold the head plate and the adjustment lever (1) of the head plate with both hands;
- Raise or lower the head plate to a suitable position;
- Release the head plate and the adjustment lever of the head plate.



The adjustment angle range of the head plate is shown in the following figure.






### Warning

- After installing the head plate make sure that the head plate is installed safely and securely to prevent the head plate from falling off and causing equipment damage or personal injury.


## 4.6 Operation of the Upper Back Plate

Press the Startup button and the Adjustment Button of the upper back plate to fold the plate up or down, or just operate and adjust the plate manually.

- 1) Firmly hold the main part of the upper back plate with one hand, hold the locking wrench (1) for upper back plate with the other hand, press the wrench button (2) with your thumb, and pull the wrench upwards;
- 2) Repeat the above operations, and loosen the locking wrench of the upper back plate on both sides;
- 3) Adjust the upper back plate up and down to a suitable position;
- 4) Pull the locking wrench downwards for upper back plate.



The adjustment angle range of the upper back plate is shown in the following figure.






### Attention

- This function is only available for W5 and W5A.

## 4.7 Operation of the Back Plate

Press the Startup button and the Adjustment Button of the back plate to fold the plate.


The adjustment angle range of the back plate is shown in the following figure.



## 4.8 Operation of the Leg Plate

### 1) Upper/lower adjustment of the leg plate

- Firmly hold the leg plate with both hands, and press the adjustment lever (1) of the leg plate.
- Raise or lower the leg plate to a suitable position.
- Release the leg plate and the adjustment lever of the leg plate.




### 2) Level adjustment of the leg plate

- Firmly hold the horizontal locking wrench (2) for the leg plate with your hand, and press the wrench button (3) with your thumb;
- Pull the horizontal locking wrench (2) for the leg plate downward;
- Adjust the leg plate left and right to a suitable position;

- d) Release the button (3) and pull up the horizontal locking wrench (2) for the leg plate to tighten the wrench.

The adjustment angle range of the leg plate is shown in the following figure.




#### Caution

- Before adjusting the leg plate make sure that there are no objects under the leg plate.

## 4.9 Lifting Operation

Press the Startup button and the Lifting Up/Down button at the same time to lift the operating table up or down.

The lifting range is shown in the following figure.




## 4.10 Forward/Backward Buckling Operation with One Key

Press the Startup button and the Forward/Backward Buckling button to adjust the operating table.


The adjustment range is shown in the following figure.






## 4.11 Operation of the Waist Elevator

Rotate the rocker (1) behind the waist elevator; this operation can raise and lower the waist elevator (2).



The adjustment range is shown in the following figure.



### Attention

- This function is only available for W3 and W3A.

## Chapter 5 Fault Handling



### Warning

- This Chapter is limited to specified users to fix simple failures. If you encounter a failure that is not included in this Chapter, or the failure cannot be eliminated after you try the fault handling methods listed below, please contact the user service organization specified by Shenzhen Comen Medical Instruments Co., Ltd. Do not repair the equipment without authorization.
- Repairs can only be done by authorized personnel of Shenzhen Comen Medical Instruments Co., Ltd. Repairs performed by unauthorized personnel may cause personal injury or equipment damage.
- Repairs must be strictly based on the technical data provided by Shenzhen Comen Medical Instruments Co., Ltd. If you need relevant technical data please contact the user service organization specified by Shenzhen Comen Medical Instruments Co., Ltd. or your local agent.

Source	Possible Cause	Countermeasures
The power plug is connected, but all light sources are off	The power switch is not turned on or the emergency stop switch is not turned on.	Turn on the emergency stop switch and the power switch.
	The wire connection is poor or faulty.	Check whether the wires in the outer cover are disconnected or fused, and take appropriate action.
	The fuse is blown.	Replace the fuse.
	Power supply is not connected.	Connect the power cord.

### 6.1 Overview

---

The equipment can be fitted with a built-in rechargeable battery. When connected to an AC power supply the battery will be charged automatically, regardless of whether the equipment is turned on or not, until it is fully charged. In the event of a sudden power failure the system will automatically use the battery to power the equipment without interrupting the equipment's status or functionality. After the AC power supply is cut off the battery indicator flashes to indicate that the battery is being used for power supply and the operation of the equipment is not affected.



#### Attention

- During storage and transportation, or when the table is not used for long periods, press the emergency stop switch of the operating table to turn off the battery output.
- If the equipment has a built-in battery the battery must be charged after each use to ensure that the battery has sufficient power reserve.
- When the battery voltage is lower than 23V and the voltage does not rise after charging it indicates that the battery has been damaged. Please replace the battery.
- When the operating table functions become intermittent or unresponsive, or the battery indicator reports a low battery indication, it indicates that the battery is too low to drive the operating table. Please charge the battery. If the behavior continues to occur after charging please contact the after-sales service.



#### Warning

- Incorrect replacement or fitting of lead-acid batteries will lead to unacceptable risks.
- Non-professional replacement of lead-acid batteries may cause risks.
- The battery electrolyte is harmful. When the battery electrolyte gets into your skin or eyes you should immediately rinse the skin or eyes with clean water and consult a physician.
- Keep the battery out of the reach of children.
- When the equipment is being powered by battery the equipment will automatically power off when the battery is low.

## 6.2 Battery Installation



### Warning

- If you need to install batteries, please contact Comen's local after-sales service.
- Only use batteries specified by the manufacturer.
- Do not remove the battery while the equipment is being used.

## 6.3 Optimization and Inspection of the Battery Performance

- Optimization of the battery performance

When the battery is used for the first time it should be subjected to at least two complete optimization cycles. A complete optimization cycle shall be uninterrupted charging until the battery is fully charged, and discharging until the equipment automatically shuts down.

When optimizing the battery, you should ensure that:

- 1) Disconnect all connections between the equipment and the patient and stop all functions.
- 2) The battery to be optimized should be fitted to the battery box.
- 3) When charging the battery, ensure that the battery is charged for more than 8 hours without interruption. The battery indicator on the side panel goes out in the shutdown state, indicating that the battery is fully charged.
- 4) Disconnect the AC power and use the battery to power the equipment until the battery is discharged and the equipment automatically shuts down.
- 5) The battery optimization is over.

- Inspection of the battery performance

The service life of the battery varies with storage, usage environment, battery discharge frequency and usage time. Even if the battery is not used its performance will gradually decrease.

The following describes the battery inspection steps:

- 1) First, determine whether the battery is damaged. When using battery power supply, if it cannot be turned on or cannot be charged (abnormal charging indicator), it means that the battery is damaged or there is no battery in the battery box.
- 2) Check whether the battery can be charged normally when the battery is connected to AC power.

- 3) Disconnect all connections between the equipment and the patient and stop all functions.
- 4) When charging the battery, ensure that the battery is charged for more than 8 hours without interruption. The battery indicator on the side panel goes out in the shutdown state, indicating that the battery is fully charged.
- 5) Disconnect the AC power and use the battery to power the equipment until the battery is discharged and the equipment automatically shuts down, and record the discharge start time and discharge end time.
- 6) The length of the discharge time reflects the battery performance.
- 7) When the discharge time drops below 50% of the initial value, the battery shall be replaced.

**Attention**

- If the battery is stored without use for long periods, in order to extend the life of the rechargeable battery, it is recommended to charge it once a month to prevent the battery from over-discharging.
- The battery life depends on the configuration and operation of the equipment. For example, frequent operation of the operating table will reduce the battery's power supply time.

## 6.4 Battery Recycling

If the battery is obviously damaged or the battery cannot be charged, it shall be replaced, and the discarded old battery shall be properly recycled. Moreover, it shall be handled according to the corresponding laws and regulations or hospital rules and regulations.

**Warning**

- Do not disassemble or short-circuit the battery or throw it into a fire, otherwise it may cause the battery to catch fire, explode, leak harmful gas or lead to other hazards.

---

## Chapter 7 Cleaning and Disinfection

---

Only use the materials and methods approved by Comen listed in this Chapter to clean or disinfect the operating table. Comen does not provide any guarantee for damage caused by using unapproved materials or methods.

Comen is not responsible for the effectiveness of the listed chemicals or methods as a means of infection control. For infection control methods, please consult your hospital's infection prevention department or epidemiologists, or refer to all local policies applicable to your hospital and country.

### 7.1 Overview

---

Please keep your equipment and accessories free of dust. After each cleaning, check the equipment carefully. If you find any signs of aging or damage, stop using it immediately. If you need to send it back to the Comen company for repairs, please clean it first. Please observe the following precautions:

- ✧ Follow the manufacturer's instructions to dilute detergents and disinfectants, or use the lowest possible concentration.
- ✧ Do not let liquid enter the housing.
- ✧ Do not let the operating table and its parts immerse in liquid for a long time.
- ✧ Do not use abrasive materials, bleaching powder, and any strong solvents (such as acetone or detergents containing acetone).



#### Warning

- **Use only the detergents and disinfectants recommended in this manual. The use of other detergents and disinfectants can cause equipment damage or lead to safety hazards.**
- **Before cleaning the operating table, you must turn off the power and disconnect the AC power supply.**
- **Do not use EtO (ethylene oxide) to disinfect the operating table.**
- **Do not let the disinfectant remain on any surface and accessories of the operating table, please wipe it off immediately with a damp cloth.**
- **Detergents cannot be mixed; otherwise dangerous gases will be generated.**
- **This Chapter only introduces the cleaning methods of reusable accessories. Disposable accessories**

cannot be cleaned and disinfected and used again to avoid cross-infection.

- In order to protect the environment, disposable accessories must be recycled or properly disposed of.
- After cleaning, if there are signs of damage or aging of the cable, the cable shall be replaced with a new one.
- The operating table and all accessories must not be autoclaved.
- Do not use any cleaning solutions other than those recommended here, which may cause permanent damage to the equipment, accessories and cables.



#### Caution

- If you accidentally pour liquid on the operating table or accessories, please contact the maintenance personnel or Comen immediately.

## 7.2 Cleaning and Disinfection of the Operating Table

The operating table shall be kept clean. It is recommended to clean the outer surface of the table frequently. In particular, in areas with severe environmental conditions such as windy or dusty, the frequency of cleaning shall be increased. In addition, in order to avoid cross-infection, it is necessary to clean the accessories regularly. Please consult or understand the hospital's regulations on equipment cleaning before cleaning.

#### ➤ Cleaning steps:

- 1) Turn off the power first and disconnect the power cord.
- 2) After using a soft cloth to absorb the proper amount of detergent, wipe the body of the operating table.
- 3) After using a soft cloth to absorb the proper amount of detergent, wipe the controller panel of the operating table.
- 4) If necessary, use a soft, dry cloth to wipe off excess detergent.
- 5) Place the equipment in a cool and ventilated environment to air dry naturally.

The disinfection operation may cause a certain degree of damage to the equipment. It is recommended that the equipment be disinfected only when necessary in the hospital maintenance plan. Before disinfecting the operating table, please clean the operating table first.

➤ Available detergents:

Cleaning/disinfection parts	Detergent	Disinfectant
Equipment housing	Ethanol (75%) Hydrogen peroxide (3%)	Isopropanol (70%)
Power cord		Glutaraldehyde solution (2%)
Fixing device of Tubes and Lines		Sodium hypochlorite disinfectant

### 7.3 Cleaning and Disinfection of Accessories

The following detergents can be used safely to clean the accessories:

General formula	Maximum concentration
Sodium hypochlorite	0.5% aqueous solution
Glutaraldehyde	2%
Hydrogen peroxide	3%



### 8.1 Maintenance inspection

---

Before the operating table is used or has been continuously used for 6 -12 months or after each maintenance or upgrade, a complete inspection of the equipment, including functional safety inspections, must be performed by trained and qualified technical maintenance personnel.

Inspection items shall include:

- 1) Check whether the operating environment and power supply of the operating table meet the requirements;
- 2) Check the operating table and accessories for mechanical damage;
- 3) Check whether the power cord is free from abrasion and has good insulation performance;
- 4) Check all equipment functions that may be used for the patient and ensure that the equipment is in a good operating condition;
- 5) Check whether the accessories used are those specified by the manufacturer;
- 6) Check whether the battery performance is good;
- 7) Check whether the wiring impedance and leakage current meet the requirements;
- 8) Confirm that the operating table has been cleaned and disinfected;
- 9) The appearance of the operating table shall be upright, uniform in color, and there shall be no obvious bumps, cracks, and sharp edges;
- 10) The texts and markings of the operating table shall be clear and accurate, and firmly stuck;
- 11) The operation of each control part shall be flexible, stable and reliable, with no jitter and no loose fasteners;
- 12) The hydraulic system of the operating table shall be leak-free.

If there is any indication of a functional impairment of the operating table is found, do not use this equipment to perform any operations on the patient. Please contact the biomedical engineer of the hospital or Comen.

All safety inspections or repairs requiring the disassembly of equipment shall be carried out by professional maintenance personnel. The use of the non-professional personnel may result in the failure of the function of the equipment or the existence of safety hazards, and may endanger personal safety.

Comen will conditionally provide circuit diagrams upon user's request to help users find appropriate and qualified technicians to repair those parts of the equipment that the Comen company classified as user-serviceable.



### Warning

- **The hospital or institution using this equipment shall establish a comprehensive maintenance plan, otherwise it may cause the failure of equipment functions and unpredictable consequences, and may endanger personal safety.**

## 8.2 Maintenance Plan

### 8.2.1 Maintenance

The following tasks can only be completed by professional maintenance personnel approved by Comen. When the following maintenance is required, please contact the maintenance personnel in time. Before testing or maintenance the operating table must be cleaned and disinfected.

Testing and maintenance items	Frequency
Carry out safety inspection according to IEC 60601-1	At least once every two years. After the equipment is dropped, or the power supply is replaced, or when required, it shall be checked.
Mechanical/hydraulic system	Carry out safety inspection and necessary maintenance once a year; The hydraulic oil shall be replaced after the first 2-years use, and then every 5 years.
Battery	See the appropriate section relating to the battery.

### 8.2.2 Service Life of Accessories

The service life of the parts and accessories of this equipment is 10 years.

---

## Appendix I Accessories

---

When this equipment is used, Comen recommends the following accessories.



### Warning

- Please use the accessory models specified by Comen. Using other accessories may damage the equipment.
- Disposable accessories can only be used once. Repeated use may result in reduced performance or cross-infection.

No.	Material Name
1	Anesthesia frame
2	Hand support
3	Leg support
4	Body/shoulder support
5	Circular clamp
6	Quadrant Clamp
7	Body strap
8	Fixing device for tubes and lines
9	Foot plate (left)
10	Foot plate (right)
11	Prone position pad

## Appendix II Product Specification

### 1) Product Type

Name	Type
CE classification	Class I
Classification in accordance with the electric shock protection type	Class I with external power supply, built-in power supply
Classification in accordance with the degree of electric shock protection	Type B
Explosion protection level	Do not provide explosion protection (general equipment), and flammable RESP agents shall not be used
EMC classification	Class A, group 1
Classification in accordance with the degree of protection against liquid ingress	This product is Class IPX4 equipment
Classification in accordance with the degree of safety when used in the case of flammable anesthetic gas mixed with air or flammable anesthetic gas mixed with oxygen or nitrous oxide	Equipment that cannot be used in the presence of flammable anesthetic gas mixed with air or flammable anesthetic gas mixed with oxygen or nitrous oxide
Classification in accordance with the operation mode	Intermittent operation, operate for 10 minutes, and rest for 20 minutes
Whether there is an applied part to protect against defibrillation discharge effect	None
Whether there is a signal output or input part	There is a signal input or output part
Permanently installed equipment or non-permanently installed equipment	Non-permanently installed equipment
Compliance with relevant standards	China: YY 0570, YY/T 1106, YY 0505, GB 9706.1, GB 4824, GB/T 4208, GB/T 191, GB/T 9969, GB/T 14710, GB/T 16886.1, GB/T 16886.5, GB/T 16886.10, YY/T 0466.1, YY/T 0316, YY/T 0664, GB/T 26572, SJ/T 11364, YY/T 1474, SJ/T 11363, Administrative Measures for the Restricted Use of Hazardous Substances in Electrical and Electronic Products, and ISTA Transport Test Procedures-3E

	CE: Medical Device Regulation (EU) 2017/745, Directive 2011/65/EU, EN ISO 14971, EN ISO 10993-1, EN ISO 10993-5, EN ISO 10993-10, ISO 780, EN 1041 ISO 15223-1, IEC 60601-1, IEC 60601-1-2, EN 60601-1-6, IEC 62366-1, IEC 60601-2-52 and IEC 60601-2-46
--	---

## 2) Physical Specifications

Name	Specification
Table surface size (excluding mattress)	Length: 2155mm $\pm$ 50mm
	Width: 540mm $\pm$ 50mm (excluding side rail)
	Height range: 680mm~980mm, error: $\pm$ 50mm
	Weight: 280kg $\pm$ 10kg

Load bearing specification	
Entire table	250kg (after the head plate and leg plate are interchanged or the longitudinal slide function is used: 190kg)

## 3) Electrical Specifications

AC power supply	
Name	Power supply environment
Power input	100 - 240V~, 50Hz/60Hz
Rated current	5.3A
Fuse	T6.3AL 250V
Power cord	5m
Built-in battery	
Voltage	24V DC (12V $\times$ 2)
Capacity	15Ah
Charging time	8h
Power supply time	About 38 operations or a week

## 4) Environmental Specifications

Working environment	
Temperature	+5°C ~ +40°C

Relative humidity	15% ~ 95%, non-condensing
Atmospheric pressure	70kPa ~ 106kPa
<b>Storage and transportation environment</b>	
Temperature	-40°C ~ +55°C
Relative humidity	10% ~ 95%, non-condensing
Atmospheric pressure	50kPa ~ 106kPa
Transport conditions	During transportation, the table must be protected from severe shock, vibration, and rain and snow.
Storage conditions	The table shall be stored in a dry and well-ventilated room with no corrosive gases.

## 5) Parameter Specifications

Intended Users	Adults, children and neonates	
Scope of application	For clinical operation	
Braking mode	Electronically controlled hydraulic brake	
Driving mode	Hydraulic electric control	
	Mechanical control	
Control mode	Hand control: Wired button control mode	
	Side control: Integrated in the column	
Data interface	Wall screen control: RS232 interface, and 1 standard interface for wall screen control	
	1 foot switch interface	
Vertical longitudinal slide		Max 320mm ±10mm
Trendelenburg / reverse Trendelenburg		Max 30° / Max 30°±5°
Lateral tilt (left / right )		Max 20° / Max 20°±5°
Fold up / fold down of the head plate		Max 60° / Max 90°±5°
Waist elevator		120mm ±10mm
Fold up / fold down of the upper back plate		Max 75° / Max 45°±5°
Fold up / fold down of the back plate		Max 80° / Max 45°±5°
Fold up / fold down of the leg plate		Max 35° / Max 90°±5°
Outward rotation of the leg plate		0° ~ 90°±5°
Lifting distance		0 ~ 300mm ±10mm
Forward/backward buckling	Forward buckling angle: 220°±5°	
	Backward buckling angle: 110°±5°	
"0" button	Trendelenburg / reverse trendelenburg, Lateral tilt, and the back plate returns to the	

	original position
The adjustment range of the table	Vertical: no more than 15mm
	Level: no more than 10mm
	Lateral: no more than 18mm
Foot control	Able to realize Trendelenburg / reverse trendelenburg, Lateral tilt and lifting function
Brake	Hydraulic brake
Movement function	The operating table is equipped with movable casters and brake; when the brake is locked, there shall be no movement at a thrust of 200N
Emergency stop	The operating table is equipped with a one-button emergency stop function
Data storage and recording	Failures data saved

## 6) Hydraulic System Specifications

Hydraulic system		
Hydraulic pump	Hydraulic working medium	Hydraulic oil (+ antioxidant and rust inhibitor) YA-N32 hydraulic oil or GB11120L-TSA32 steam turbine oil
	Rated working pressure	21Mpa
	Operating working pressure range	10Mpa±10%
Hydraulic Leakage	When hydraulic leakage lasts for 5 hours, the drop distance of the table should not exceed 8mm.	
Connector and installation method	Hose connector	Withholding type hose connector (with good pulling resistance and sealing performance)
	Installation method	Threaded (tubular) mounting connection

## 7) Warning

Low battery shutdown delay time	5 minutes after the first low battery warning
---------------------------------	---

---

## Appendix III EMC

---



### Attention

- The W5, W5A, W3 and W3A operating table meets the requirements of electromagnetic compatibility in the IEC 60601-1-2 and IEC 60601-2-46 standards.
- Users shall install and use the electromagnetic compatibility information provided in the accompanying file.
- Portable and mobile RF communication equipment may affect the performance of the W5, W5A, W3 and W3A operating table. Avoid strong electromagnetic interference during use, such as near mobile phones, microwave ovens, etc.
- See the attachments below in the manual for the guidance and manufacturer's declaration.



### Warning

- The W5, W5A, W3 and W3A operating tables shall not be used close to other equipment. If they must be used close to each other, all devices shall be observed to verify that they can operate normally in the configuration used.
- Class A equipment is intended for use in industrial environments. Due to conducted emissions and radiated emissions in the W5, W5A, W3 and W3A operating table, it may be potentially difficult to ensure electromagnetic compatibility in other environments.
- In addition to the cables sold by the manufacturer of the W5, W5A, W3 and W3A operating table as spare parts for internal components, the use of accessories and cables outside the regulations may result in increased emission or reduced immunity of the W5, W5A, W3 and W3A operating table.
- Even if other equipment meet the emission requirements of corresponding national standards, the equipment or system may still be interfered by other equipment.



Attachments:

<b>Cable information</b>				
<b>No.</b>	<b>Name</b>	<b>Cable length (m)</b>	<b>Shield Radio Signal or not</b>	<b>Remark</b>
1	AC power cord	5	No	/
2	Hand control cable	3.5	No	/

Table 1


<b>Guidance and manufacturer's declaration - electromagnetic emission</b>		
W5, W5A, W3 and W3A operating table is intended for use in the electromagnetic environment specified below. The customer or the user of W5, W5A, W3 and W3A operating table should assure that it is used in such an environment.		
<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF conduction CISPR 11	Group 1	The W5, W5A, W3 and W3A operating tables use RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Compliant	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Compliant	

Table 2

<b>Guidance and manufacturer's declaration - electromagnetic immunity</b>			
The W5, W5A, W3 and W3A operating tables are intended for use in the electromagnetic environment specified below. The customer or user of W5, W5A, W3 and W3A operating table should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	$\pm 8$ kV contact $\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines 100 kHz repetition frequency	$\pm 2$ kV for power supply lines $\pm 1$ kV for input/output lines 100 kHz repetition frequency	The mains power quality should meet the requirements of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 0.5$ kV, $\pm 1$ kV line(s) to lines $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV line(s) to earth	$\pm 0.5$ kV, $\pm 1$ kV line(s) to lines $\pm 0.5$ kV, $\pm 1$ kV, $\pm 2$ kV line(s) to earth	The mains power quality should meet the requirements of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$>95\% U_T$ ; 0.5 cycle At $0^\circ, 45^\circ, 90^\circ, 135^\circ, 180^\circ, 225^\circ, 270^\circ$ and $315^\circ$ $>95\% U_T$ ; 1 cycle At $0^\circ$ $30\% U_T$ ; 25/30 cycles Single phase: at $0^\circ$ $>95\% U_T$ ; 250/300 cycles	$>95\% U_T$ ; 0.5 cycle At $0^\circ, 45^\circ, 90^\circ, 135^\circ, 180^\circ, 225^\circ, 270^\circ$ and $315^\circ$ $>95\% U_T$ 1 cycle At $0^\circ$ $30\% U_T$ ; 25/30 cycles Single phase: at $0^\circ$ $>95\% U_T$ ; 250/300 cycles	The mains power quality should meet the requirements of a typical commercial or hospital environment. If the user of the W5, W5A, W3 and W3A operating table requires continued operation during power interruptions, it is recommended that the W5, W5A, W3 and W3A operating table be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field IEC 61000-4-8	30 A/m, 50/60 Hz	30 A/m	Power frequency magnetic fields should have the magnetic characteristic of a typical location in a typical commercial or hospital

			environment.
NOTE: $U_T$ is the a.c. mains voltage prior to the application of the test level.			

Table 3

Guidance and manufacturer declaration - electromagnetic immunity			
The W5, W5A, W3 and W3A operating tables are intended for use in the electromagnetic environment specified below. The customer or the user of the W5, W5A, W3 and W3A operating table should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V RMS 0.15 MHz to 80 MHz 6 V in ISM and between 0.15 MHz and 80 MHz	3 V RMS 0.15 MHz to 80 MHz 6 V in ISM and between 0.15 MHz and 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the W5, W5A, W3 and W3A operating table, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz	3V/m	<p>Recommended separation distance</p> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P} \quad 150 \text{ KHz to } 80 \text{ MHz}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.7 \text{ GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup></p> <p>Interference may occur near the equipment marked with the following symbol:</p> 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection from structures, objects and people.

- a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which W5, W5A, W3 and W3A operating table is used exceeds the applicable RF compliance level above, W5, W5A, W3 and W3A operating table should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating W5, W5A, W3 and W3A operating table.
- b Over the frequency range 0.15 MHz to 80 MHz, field strengths should be less than 3 V/m.

Table 4

Recommended separation distances between portable and mobile RF communications equipment and W5, W5A, W3 and W3A Operating Table			
The W5, W5A, W3 and W3A operating table is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of W5, W5A, W3 and W3A operating tables can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the W5, W5A, W3 and W3A operating table, as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	0.15 MHz to 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.7 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range is applied.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by the			

absorption and reflection from structures, objects and people.

Table 5

Declaration - IMMUNITY to proximity fields from RF wireless communications equipment						
W5, W5A, W3 and W3A operating table is intended for use in an electromagnetic environment in which RF wireless communications equipment are controlled.						
Immunity test	IEC60601 test level				Compliance level	Electromagnetic environment - guidance
	Test frequency	Modulation	Maximum power	Immunity level		
Radiated RF IEC 61000-4-3	385 MHz	**Pulse Modulation: 18Hz	1.8W	27 V/m	27 V/m	
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	28 V/m	
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m	
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18Hz	2 W	28 V/m	28 V/m	
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m	
	2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m	
	5240 MHz 5500 MHz 5785 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m	
	Note * - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.					
Note** - The carrier shall be modulated using a 50 % duty cycle square wave signal.						



## Appendix IV Toxic and Hazardous Substances or Elements in the Product

Part name		Lead Pb	Mercury Hg	Cadmium Cd	Hexavalent chromium Cr (VI)	Polybrominated biphenyls PBB	Polybrominated diphenyl ethers PBDE
Equipment housing	Metal element	×	0	0	×	0	0
	Table surface	0	0	0	0	0	0
	FRP cover	0	0	0	×	0	0
	Label	0	0	0	0	0	0
Main unit	column	0	0	0	×	0	0
	Main unit hardware	0	0	0	×	0	0
	Cylinder	0	0	0	×	0	0
	Pipeline	0	0	0	0	0	0
	Solenoid valve	0	0	0	0	0	0
	Fuel pump	0	0	0	0	0	0
	Internal cable	0	0	0	0	0	0
	PCBA	×	0	0	0	0	0
Packaging	Packaging material	×	×	0	0	×	×
Universal	Connector	0	0	0	×	0	0
	Power cord	0	0	0	0	0	0
Battery	Lead-acid battery	×	×	×	×	×	×
Accessories	Metal element	×	0	0	×	0	0
	Sheet	0	0	0	0	0	0
	Pad	0	0	0	0	0	0
Remark	<p>0: It indicates that the content of this toxic and hazardous substance in all homogeneous materials of this part is below the limit requirement stipulated in SJ/T 11363-2006.</p> <p>×: It indicates that the content of this toxic and hazardous substance in at least one homogeneous material of this part exceeds the limit requirement stipulated in SJ/T 11363-2006.</p>						

---

## Appendix V Consideration for Environmentally Conscious Design

---

### 1. Instructions for Minimizing Environmental Impact during Normal Use

This part is compiled based on the requirements of *Clause 4 Protection of Environment, 4.5.2 Instructions for minimizing environmental impact during normal use* of IEC 60601-1-9.

According to the requirements of this clause, the manufacturer shall provide instructions for minimizing the environmental impact of the ME equipment during normal use in the accompanying documents.

The instructions cover the following items (Table 1).

**Table 1 The requirements of Clause 4.5.2 and Instructions provided by manufacturer**

The requirements of Clause 4.5.2	Instructions provided by manufacturer
1) Instructions on how to install the ME EQUIPMENT in order to minimize the ENVIRONMENTAL IMPACT during its EXPECTED SERVICE LIFE;	Try to keep the integrity of the non-disposable packing materials and put away the packing materials for future use or put into the specified location where complying with the rules and regulations of the local government and the hospital. Avoid overusing the cleaning reagents and other substances. For the reusable accessories, clean them with specified reagent and put away, and for the disposable ones, deal with them together and put into the specified location where complying with the rules and regulations of the local government and the hospital. If not specified, please follow the rules and regulations of the local government and the hospital.
2) Instructions on how to use and maintain the ME EQUIPMENT in order to minimize the ENVIRONMENTAL IMPACT during its EXPECTED SERVICE LIFE;	Use the specified accessories and cleaning and disinfection reagent to avoid harm to the machine and accessories and reduction of the service life. Use the medical device strictly according to the instruction manual. And for maintaining the medical device, always dilute according to the manufacturer's instructions or use lowest possible concentration. Never use bleach. Do not mix disinfecting solutions (such as bleach and ammonia) as this may result in hazardous or poisonous gases or liquids. When there is a need to maintain, please follow the Instruction for Use or follow the rules and regulations of the Hospital.
3) Consumption during NORMAL USE (e.g. energy, consumable materials/parts,	During normal use of this device, it will consume electricity (alternate current and direct current-battery).



disposables, water, gasses, chemicals/reagents etc.);	The disposables are also consumed and shall be disposed according to the rules of the local government and the hospital. For cleaning or disinfection for the cables and machine, the water and ethanol or isopropanol will be used and the waste liquid shall be thrown according to the rules of the local government and the hospital.
4) Emissions during NORMAL USE (e.g. WASTE water, WASTE consumable materials, acoustic energy, heat, gasses, vapours, particulates, HAZARDOUS SUBSTANCES and other WASTE);	During normal use, it is expected there will be some consumption of the medical device. To avoid unnecessary consumption such as acoustic energy, heat, gases, hazardous substances, etc., it's recommended that on the premise of normal operation, turn off the unused functions to reduce the unnecessary heat emission and electricity consumption.
5) Information on the location within the ME EQUIPMENT of HAZARDOUS SUBSTANCES, radioactive sources and induced radioactive materials.	The battery is located inside the table column. Capacitors may contain stored energy or may pose other hazards, assembled on the PCB boards within the device.

## 2. Information for End of Life Management

This part is compiled based on *Clause 4 Protection of Environment, 4.5.3 Information for end of life management* of IEC 60601-1-9.

According to the requirements of this clause, the manufacturer shall provide the responsible organization with information for the proper disposal of the ME equipment at *End of Life* (EOL). And the manufacturer shall make available information to waste treatment facilities necessary for the environmentally responsible management of end of life ME equipment.

The information shall contain the following items (Table 2).

**Table 2 The requirements of Clause 4.5.3 and Instructions provided by manufacturer**

The requirements of Clause 4.5.3	Instructions provided by manufacturer
1) The location of components and parts within the ME equipment that contain stored energy or pose other hazards that can result in an unacceptable risk to disassemblers or others and methods for controlling such risks.	The battery is located inside the table column. Capacitors may contain stored energy or may pose other hazards, assembled on the PCB boards within the device.
2) The identity and location of hazardous substances requiring special handling and	The battery is located inside the table column. Capacitors may contain stored energy or may pose other hazards, assembled on the PCB boards within the

treatment	device.
<p>3) Disassembly instructions sufficient for the safe removal of these hazardous substances including radioactive sources and induced radioactive materials within the ME equipment.</p>	<p>For other hazards that may result in unacceptable risk, the main concern is the handling with battery: Risk of fire, explosion, or burns. Do not crush, puncture, disassemble or short circuit the battery. Do not dispose of the battery in fire or water. Do not place the battery in an environment whose temperature is above 60°C (140°F). Store the battery in the -20°C (-4°F) to 60°C (140°F) environment. Use the specified charger only. Read instructions for use. Maximum Recommended Ambient temperature is 45°C (125°F).</p> <p>Dispose of used batteries promptly and in an environmentally-responsible manner. Do not dispose of the battery in normal waste containers. Consult your hospital administrator to find out about local arrangements.</p> <p>As for disposing of the medical device, to avoid contaminating or infecting personnel, the environment or other equipment, make sure you disinfect and decontaminate the medical device appropriately before disposing of it in accordance with your country's laws for equipment containing electrical and electronic parts. For disposal of parts and accessories such as thermometers, where not otherwise specified, follow local regulations regarding disposal of hospital waste.</p>